

# 2 Designs. 3 Gauges. 12 Models. There's a SurgiLance™ for Everybody!



✓	EASY TO USE	✓
✓	MINIMIZES PAIN	✓
✓	SAFE, SINGLE-USE DESIGN	✓
✓	PROVEN RELIABILITY	✓
✓	21 GAUGE NEEDLES	✓
	28 GAUGE NEEDLES	✓
✓	18 GAUGE BLADES	
	ECO-FRIENDLY DESIGN	✓
	EASIER FOR SMALL HANDS	✓

With **more than a billion devices sold since 2000**, there's a reason why the SurgiLance™ brand from MediPurpose™ has remained a **top choice for capillary blood collection fingerstick users**: They love SurgiLance's **ease of use, comfort, safety and reliability!**

And because it's a medical product from MediPurpose, customers know they are getting **excellent value** from a company known for its commitments to **quality and customer service.**

[medipurpose.com/surgilance](http://medipurpose.com/surgilance)

**MediPurpose™**

# 2 Designs. 3 Gauges. 12 Models. There's a SurgiLance™ for Everybody!



MODEL	GAUGE†	DEPTH	FLOW	USE
SLN100	G21	1.0 mm	low 5-10 µL	glucose test geriatrics
SLN200	G21	1.8 mm	low 10-20 µL	glucose test
SLN240	G21	2.2 mm	medium 20-40 µL	glucose test hematocrit tube
SLN300	G21	2.8 mm	high 40-60 µL	glucose test hematocrit tube
SLB200	G18	1.8 mm	high 70-100 µL	microcollection tube
SLB250	G18	2.3 mm	high 100-150 µL	microcollection tube

†Standard Wire Gauge



MODEL	GAUGE†	DEPTH	FLOW	USE
SLL140	G28	1.4 mm	low 2-5 µL	glucose test
SLL180	G28	1.8 mm	low 5-10 µL	glucose test
SLL220	G28	2.2 mm	low 10-20 µL	glucose test
SLL21G180	G21	1.8 mm	low 10-20 µL	glucose test
SLL21G220	G21	2.2 mm	medium 20-40 µL	glucose test hematocrit tube
SLL21G280	G21	2.8 mm	high 40-60 µL	glucose test hematocrit tube

[medipurpose.com/surgilance](http://medipurpose.com/surgilance)

