

# 3 Designs. 3 Gauges. 13 Models. There's a SurgiLance™ for Everybody!



EASY TO USE	✓	✓	✓
MINIMIZES PAIN	✓	✓	✓
SAFETY LANCET	✓		✓
PROVEN RELIABILITY	✓	✓	✓
21 GAUGE NEEDLES	✓		✓
28 GAUGE NEEDLES	✓	✓	✓
18 GAUGE BLADES	✓		
ECO-FRIENDLY DESIGN			✓
EASIER FOR SMALL HANDS			✓

With **more than a billion devices sold since 2000**, there's a reason why the SurgiLance™ brand from MediPurpose™ has remained a **top choice for capillary blood collection fingerstick users**: They love SurgiLance's **ease of use, comfort, safety and reliability!**

And because it's a medical product from MediPurpose, customers know they are getting **excellent value** from a company known for its commitments to **quality and customer service.**

[medipurpose.com/surgilance](https://www.medipurpose.com/surgilance)

**MediPurpose™**

# 3 Designs. 3 Gauges. 13 Models. There's a SurgiLance™ for Everybody!

MODEL	GAUGE†	DEPTH	FLOW	USE
SLN170	G28	1.7 mm	low 5–10 µL	glucose test geriatrics
SLN100	G21	1.0 mm	low 5–10 µL	glucose test geriatrics
SLN200	G21	1.8 mm	low 10–20 µL	glucose test
SLN240	G21	2.2 mm	medium 20–40 µL	glucose test hematocrit tube
SLN300	G21	2.8 mm	high 40–60 µL	glucose test hematocrit tube
SLB200	G18	1.8 mm	high 70–100 µL	microcollection tube
SLB250	G18	2.3 mm	high 100–150 µL	microcollection tube



† Standard Wire Gauge

MODEL	GAUGE†	FLOW	USE††
SLR28G	G28	low 2–5 µL	glucose test



†† Fits most lancing devices for capillary blood sampling

MODEL	GAUGE†	DEPTH	FLOW	USE
SLL28G140	G28	1.4 mm	low 2–5 µL	glucose test
SLL28G180	G28	1.8 mm	low 5–10 µL	glucose test
SLL28G220	G28	2.2 mm	low 10–20 µL	glucose test
SLL21G180	G21	1.8 mm	low 10–20 µL	glucose test
SLL21G220	G21	2.2 mm	medium 20–40 µL	glucose test hematocrit tube
SLL21G280	G21	2.8 mm	high 40–60 µL	glucose test hematocrit tube



[medipurpose.com/surgilance](https://www.medipurpose.com/surgilance)

**MediPurpose™**