



Multipurpose capillary blood sampling devices for all testing requirements.

Safe, easy and effective solutions for fingersticks!



CAT. No. AND TYPE	COLOR	DEPTH OF PENETRATION	GAUGE [†] SWG	BLOOD FLOW (APPROXIMATE)	USE
SLN100 Needle	Yellow	1.0 mm	G21	Low Flow 5–10 µL	Glucose Test Geriatrics
SLN170 Needle	Lime	1.7 mm	G28	Low Flow 5–10 µL	Glucose Test
SLN200 Needle	Gray	1.8 mm	G21	Low Flow 10–20 µL	Glucose Test
SLN240 Needle	Orange	2.2 mm	G21	Medium Flow 20–40 µL	Glucose Test Hematocrit Tube
SLN300 Needle	Pink	2.8 mm	G21	High Flow 40–60 µL	Glucose Test Hematocrit Tube
SLB200 Blade	Green	1.8 mm	G18	High Flow 70–100 µL	Microcollection Tube
SLB250 Blade	Blue	2.3 mm	G18	High Flow 100–150 µL	Microcollection Tube

[†]SWG: Standard Wire Gauge at the point of penetration.

Americas | United States Asia | Singapore

2862 Buford Hwy
Suite 105
Duluth, Georgia 30096, USA
Tel: +1-404-855-3655
Email: sales@medipurpose.com

10 Anson Road
#12-08 International Plaza
Singapore 079903
Email: sales@medipurpose.com

Europe | United Kingdom

3 College Gardens
New Malden, Surrey KT3 6NT
England, United Kingdom
Email: sales@medipurpose.com

MediPurpose™
medipurpose.com/surgilance