



Competitive Safety Lancet Cross Reference Guide*

	Competitive Product					SurgiLance® Equivalent			
	Brand	Product #	Color	Needle/ Blade	Specifications	Product #	Color	Needle/ Blade	Specifications
ACP MEDICAL (FORMERLY VITAL CARE) 	Auto Junior Auto Soft Auto Normal Auto Extra	FTM065 FTM066 FTM067 FTM069	Pink Blue Beige Green	Needle Needle Needle Needle	26 gauge x 1.8 mm depth 21 gauge x 1.8 mm depth 21 gauge x 2.4 mm depth 21 gauge x 3.0 mm depth	SLN100 SLN300 SLN300 SLN300	Yellow Pink Pink Pink	Needle Needle Needle Needle	26 gauge x 1.0 mm depth 21 gauge x 2.8 mm depth 21 gauge x 2.8 mm depth 21 gauge x 2.8 mm depth
ACTI-LANCE 	Acti-Lance Lite Acti-Lance Universal Acti-Lance Special	7141 7142 7143	Purple Blue Gold	Needle Needle Blade	28 gauge x 1.5 mm depth 23 gauge x 1.8 mm depth 17 gauge x 2.0 mm depth	SLN170 SLN200 SLB200	Lime Grey Blue	Needle Needle Blade	28 gauge x 1.7 mm depth 23 gauge x 1.8 mm depth 18 gauge x 2.3 mm depth
ARKRAY/HYPOGUARD 	Haemolance* Haemolance* Haemolance Plus Haemolance Plus Haemolance Plus Haemolance Plus Haemolance Plus Haemolance Plus Haemolance Plus Haemolance Plus *	990300 990500 990200 990900 990700 990800 990400 990600	Blue Green Light Blue Blue Green Yellow Pink Purple	Needle Needle Needle Needle Needle Needle Blade Blade	21 gauge x 1.8 mm depth 25 gauge x 1.8 mm depth 28 gauge x 1.6 mm depth 25 gauge x 1.4 mm depth 21 gauge x 1.8 mm depth 18 gauge x 1.9 mm depth 1.2 mm depth 1.6 mm depth	SLN300 SLN100 SLN170 SLN100 SLN300 SLB250 SLB200 SLB250	Pink Yellow Lime Yellow Pink Blue Green Blue	Needle Needle Needle Needle Needle Blade Blade Blade	21 gauge x 2.8 mm depth 26 gauge x 1.0 mm depth 28 gauge x 1.7 mm depth 26 gauge x 1.0 mm depth 21 gauge x 2.8 mm depth 18 gauge x 2.3 mm depth 19 gauge x 1.8 mm depth 18 gauge x 2.3 mm depth
BD 	Genie Genie Genie Genie Genie QuikHeel QuikHeel	366579 366583 366580 366581 366582 368100 368101	Purple Orange Pink Green Blue Purple Teal	Needle Needle Blade Blade Blade Blade Blade	28 gauge x 1.25 mm depth 23 gauge x 2.25 mm depth 1.0 mm x 1.50 mm depth 1.5 mm x 1.50 mm depth 2.0 mm x 1.50 mm depth 1.75 mm x 0.85 mm depth 2.5 mm x 1.00 mm depth	SLN170 SLN200 SLB200 SLB200 SLB250 SLB200 SLB250	Lime Grey Green Green Blue Green Blue	Needle Needle Blade Blade Blade Blade Blade	28 gauge x 1.7 mm depth 23 gauge x 1.8 mm depth 19 gauge x 1.8 mm depth 19 gauge x 1.8 mm depth 18 gauge x 2.3 mm depth 19 gauge x 1.8 mm depth 18 gauge x 2.3 mm depth
CH DIAGNOSTICS 	Novalance	4225	Blue	Needle	25 gauge x 1.8 mm depth	SLN100	Yellow	Needle	26 gauge x 1.0 mm depth
FISHER SCIENTIFIC 	Fisherbrand Low Fisherbrand Med Fisherbrand High	02-675-161 02-675-160 02-675-165	Blue Green Yellow w	Needle Needle Needle	25 gauge x 1.8 mm depth 21 gauge x 1.8 mm depth 18 gauge x 1.8 mm depth	SLN100 SLN300 SLB250	Yellow Pink Blue	Needle Needle Blade	26 gauge x 1.0 mm depth 21 gauge x 2.8 mm depth 18 gauge x 2.3 mm depth
HAWAII MEDICAL 	NeatNick NeatNick	1030083 1030085	Lime Green Mango Orange	Blade Blade	1.0 mm depth 1.2 mm depth	SLB200 SLB200	Green Green	Blade Blade	0.92 mm x 1.8 mm depth 0.92 mm x 1.8 mm depth

*NOTE: Items are cross-referenced by needle size first, then by depth. Cross-reference specifications were obtained via each manufacturer's website in October 2007, however, data is subject to change without notice.



Competitive Safety Lancet Cross Reference Guide*

	Competitive Product					SurgiLance® Equivalent			
	Brand	Product #	Color	Needle/ Blade	Specifications	Product #	Color	Needle/ Blade	Specifications
HELENA LABORATORIES 	Little Nicky Nicky	LNK50 NK50	White Blue	Blade Blade	0.85 mm x 1.75 mm depth 1.00 mm x 2.50 mm depth	SLB200 SLB250	Green Blue	Blade Blade	0.92 mm x 1.8 mm depth 1.15 mm x 2.3 mm depth
INVACARE 	Low Flow Normal Flow	ISG171232 ISG171231	Blue Green	Needle Needle	25 gauge 23 gauge	SLN100 SLN200	Yellow Grey	Needle Needle	26 gauge x 1.0 mm depth 23 gauge x 1.8 mm depth
ITC 	Tenderlett Tenderlett Tenderlett Tenderfoot Tenderfoot Tenderfoot Tenderfoot	TL100 TLJ100 TLT100 TFM50 TFP50 TF50 TFT50	Red (adult) Blue (junior) Pink (toddler) Blue White Blue/Pink Pink	Blade Blade Blade Blade Blade Blade Blade	1.75 mm x 0.94 mm depth 1.25 mm x 0.67 mm depth 0.85 mm x 0.46 mm depth 0.65 mm x 1.40 mm depth 0.85 mm x 1.75 mm depth 1.00 mm x 2.50 mm depth 2.00 mm x 3.00 mm depth	SLB200 SLB200 SLB200 SLB200 SLB200 SLB250 SLB250	Green Green Green Green Green Blue Blue	Blade Blade Blade Blade Blade Blade Blade	0.92 mm x 1.8 mm depth 0.92 mm x 1.8 mm depth 0.92 mm x 1.8 mm depth 0.92 mm x 1.8 mm depth 0.92 mm x 1.8 mm depth 1.15 mm x 2.3 mm depth 1.15 mm x 2.3 mm depth
LIFESCAN 	Tenderfoot One-Touch SureSoft Gentle One-Touch SureSoft Regular One-Touch SureSoft Dual	021-140-01 021-139-01 021-141-01	Light Blue Dark Blue Yellow	Needle Needle Needle	28 gauge x 1.8 mm depth 21 gauge x 2.0 mm depth 18 gauge x 1.8 mm depth	SLN170 SLN300 SLB250	Lime Pink Blue	Needle Needle Blade	28 gauge x 1.7 mm depth 21 gauge x 2.8 mm depth 18 gauge x 2.3 mm depth
LIFESCAN (FORMERLY OWN MUMFORD) 	Unistik 2 Neonatal Unistik 2 Regular Unistik 2 Deep	020-196 010-636 010-404	Yellow Light Blue Dark Blue	Blade Needle Needle	18 gauge x 1.8 mm depth 21 gauge x 2.4 mm depth 21 gauge x 3.0 mm depth	SLB200 SLN300 SLN300	Green Pink Pink	Blade Needle Needle	19 gauge x 1.8 mm depth 21 gauge x 2.8 mm depth 21 gauge x 2.8 mm depth
MEDIPURPOSE 	SurgiLance Safety Lancet SurgiLance Safety Lancet SurgiLance Safety Lancet SurgiLance Safety Lancet SurgiLance Safety Lancet SurgiLance Safety Lancet SurgiLance Safety Lancet	SLN100 SLN170 SLN200 SLN240 SLN300 SLB200 SLB250	Yellow Lime Grey Orange Pink Green Blue	Needle Needle Needle Needle Needle Blade Blade	26 gauge x 1.0 mm depth 28 gauge x 1.7 mm depth 23 gauge x 1.8 mm depth 22 gauge x 2.2 mm depth 21 gauge x 2.8 mm depth 19 gauge x 1.8 mm depth 18 gauge x 2.3 mm depth				
OWEN MUMFORD 	Unistik 3 Comfort Unistik 3 Normal Unistik 3 Extra Unistik 3 Dual	AT104 2 AT100 2 AT101 2 AT106	Lavender Yellow Orange Blue	Needle Needle Needle Blade	28 gauge x 1.8 mm depth 23 gauge x 1.8 mm depth 21 gauge x 2.0 mm depth 18 gauge x 1.8 mm depth	SLN170 SLN200 SLN300 SLB250	Lime Grey Pink Blue	Needle Needle Needle Blade	28 gauge x 1.7 mm depth 23 gauge x 1.8 mm depth 21 gauge x 2.8 mm depth 18 gauge x 2.3 mm depth
*NOTE: Items are cross-referenced by needle size first, then by depth. Cross-reference specifications were obtained via each manufacturer's website in October 2007, however, data is subject to change without notice.									

Competitive Safety Lancet Cross Reference Guide*

	Competitive Product					SurgiLance® Equivalent			
	Brand	Product #	Color	Needle/ Blade	Specifications	Product #	Color	Needle/ Blade	Specifications
PALCO LABS 	Super Soft Soft Normal Ultra	ELY100 ELG100 ELO100 ELR100	Yellow Green Orange Blue	Needle Needle Needle Blade	26 gauge x 1.8 mm depth 21 gauge x 1.8 mm depth 21 gauge x 2.2 mm depth 21 gauge x 2.2 mm depth	SLN100 SLN300 SLN300 SLB200	Yellow Pink Pink Green	Needle Needle Needle Blade	26 gauge x 1.0 mm depth 21 gauge x 2.8 mm depth 21 gauge x 2.8 mm depth 19 gauge x 1.8 mm depth
ROCHE 	Safe-T-Pro	951	Blue/White	Needle	1.8 mm depth	SLN200	Grey	Needle	23 gauge x 1.8 mm depth
SMITHS MEDICAL (FORMERLY FUTURA MEDICAL) 	Micro Flow Low Flow Normal Flow High Flow Maximum Flow Pediatric Blade	1025 1020 1014 1008 1005 1028	Light Blue Blue Green Yellow Purple Pink	Needle Needle Needle Needle Blade Blade	28 gauge x 1.8 mm depth 25 gauge x 1.8 mm depth 21 gauge x 1.8 mm depth 18 gauge x 1.8 mm depth	SLN170 SLN100 SLN300 SLB250 SLB250 SLB200	Lime Yellow Pink Blue Blue Green	Needle Needle Needle Blade Blade Blade	28 gauge x 1.7 mm depth 26 gauge x 1.0 mm depth 21 gauge x 2.8 mm depth 18 gauge x 2.3 mm depth 18 gauge x 2.3 mm depth 19 gauge x 1.8 mm depth
STAT MEDICAL 	Pediatric/Tender Comfort Normal Neonatal/Heelstick PediLet PediLet Preemie	SAP 100 SAC 100 SAN 100 SAH 100 SPS 50 SPP 50	Light Blue Grey Orange Dark Blue Blue/Pink Pink/White	Needle Needle Needle Blade Blade Blade	26 gauge x 1.8 mm depth 21 gauge x 1.8 mm depth 21 gauge x 2.4 mm depth 18 gauge x 1.2 mm depth 1.00 mm x 2.50 mm depth 0.85 mm x 1.75 mm depth	SLN100 SLN300 SLN300 SLB250 SLB250 SLB200	Yellow Pink Pink Blue Blue Green	Needle Needle Needle Blade Blade Blade	26 gauge x 1.0 mm depth 21 gauge x 2.8 mm depth 21 gauge x 2.8 mm depth 18 gauge x 2.3 mm depth 18 gauge x 2.3 mm depth 19 gauge x 1.8 mm depth
TYCO HEALTHCARE 	Monolettor	88816020091	Blue/White	Needle	21 gauge x 2.4 mm depth	SLN300	Pink	Needle	21 gauge x 2.8 mm depth

*NOTE: Items are cross-referenced by needle size first, then by depth. Cross-reference specifications were obtained via each manufacturer's website in October 2007, however, data is subject to change without notice.