



THIS  
IS A  
**babyLance™**  
BABY

babyLance™ Safety Heelstick Product Guide

**Designed for neonatal nurses...with neonatal nurses!**

babyLance™ safety heelsticks have earned worldwide popularity for their **smooth and accurate incisions**, **ease of use** and **safety**.





## Designed for neonatal nurses...with neonatal nurses!

babyLance™ safety heelsticks have earned worldwide popularity for their **smooth and accurate incisions, ease of use** and **safety**.

### Solutions for Micropreemies to Toddlers

After continuing to listen to the needs of neonatal caregivers, we have expanded the babyLance™ family to include **new Micropreemie and Toddler devices**.

### Ultra-Smooth Incisions Minimize Trauma

babyLance's proprietary internal spring activates its blade's swift pendulum action—delivering an **ultra-smooth incision** that complies with CLSI LA4-A5 guidelines to **reduce trauma** to an infant's delicate subcutaneous tissues.





### Safety

**Single-use and auto-retraction** design practically eliminates cross-contamination and accidental needlestick injuries—before or after use.

### Ease of Use

Featuring the industry's most **effortless activation trigger**, babyLance's unique housing provides visual guides for **accuracy** and dimples for a **secure grip**.

## A Preferred Safety Heelstick for Every Neonatal Heel Incision Procedure

	babyLance™ Product	Incision Depth	Indications <sup>†</sup>
	Micropreemie (BLM)	0.60 mm	Less Than 1000 Grams <sup>†</sup>
	Premie (BLP)	0.85 mm	1000–2500 Grams Low Birth Weight <sup>†</sup>
	Newborn (BLN)	1.0 mm	2500 Grams–9 Kilograms Birth–Six Months <sup>†</sup>
	Toddler (BLT)	1.4 mm	More Than 9 Kilograms Six Months–Two Years <sup>†</sup>

<sup>†</sup> These are suggested indications; the medical staff's clinical judgement should determine which babyLance™ model to use.

#### Americas

3883 Rogers Bridge Road NW  
Suite 501  
Duluth, GA 30097  
Tel: +1 770 448 9493

#### Asia

15 Hoe Chiang Road #12-02  
Tower Fifteen  
Singapore 089316  
Tel: +65 63451588

#### Europe

Tel: +44 330 122 2630

**MediPurpose™**

medipurpose.com/babyLance